

MULTIPLES OF ELEMENT WEIGHTS

C	12.01	H ₈	8.064	H ₅₀	50.40	Cl ₄	141.81
C ₂	24.02	H ₉	9.072	H ₅₁	51.41	Cl ₅	177.27
C ₃	36.03	H ₁₀	10.08	H ₅₂	52.42	Br	79.90
C ₄	48.04	H ₁₁	11.09	H ₅₃	53.42	Br ₂	159.81
C ₅	60.05	H ₁₂	12.10	H ₅₄	54.43	Br ₃	239.71
C ₆	72.06	H ₁₃	13.10	H ₅₅	55.44	Br ₄	319.62
C ₇	84.07	H ₁₄	14.11	H ₅₆	56.45	I	126.90
C ₈	96.08	H ₁₅	15.12	H ₅₇	57.46	I ₂	253.80
C ₉	108.09	H ₁₆	16.13	H ₅₈	58.46	I ₃	380.70
C ₁₀	120.10	H ₁₇	17.14	H ₅₉	59.47	OCH ₃	31.03
C ₁₁	132.11	H ₁₈	18.14	H ₆₀	60.48	(OCH ₃) ₂	62.07
C ₁₂	144.12	H ₁₉	19.15	H ₆₁	61.49	(OCH ₃) ₃	93.10
C ₁₃	156.13	H ₂₀	20.16	H ₆₂	62.50	(OCH ₃) ₄	124.14
C ₁₄	168.14	H ₂₁	21.17	H ₆₃	63.50	(OCH ₃) ₅	155.17
C ₁₅	180.15	H ₂₂	22.18	H ₆₄	64.51	(OCH ₃) ₆	186.20
C ₁₆	192.16	H ₂₃	23.18	H ₆₅	65.52	(OCH ₃) ₇	217.24
C ₁₇	204.17	H ₂₄	24.19	O	16.00	(OCH ₃) ₈	248.27
C ₁₈	216.18	H ₂₅	25.20	O ₂	32.00	OC ₂ H ₅	45.06
C ₁₉	228.19	H ₂₆	26.21	O ₃	48.00	(OC ₂ H ₅) ₂	90.12
C ₂₀	240.20	H ₂₇	27.22	O ₄	64.00	(OC ₂ H ₅) ₃	135.18
C ₂₁	252.21	H ₂₈	28.22	O ₅	80.00	(OC ₂ H ₅) ₄	180.24
C ₂₂	264.22	H ₂₉	29.23	O ₆	96.00	(OC ₂ H ₅) ₅	225.30
C ₂₃	276.23	H ₃₀	30.24	O ₇	112.00	OCOCH ₃	59.04
C ₂₄	288.24	H ₃₁	31.25	O ₈	128.00	(OCOCH ₃) ₂	118.09
C ₂₅	300.25	H ₃₂	32.26	O ₉	144.00	(OCOCH ₃) ₃	177.13
C ₂₆	312.26	H ₃₃	33.26	O ₁₀	160.00	(OCOCH ₃) ₄	236.18
C ₂₇	324.27	H ₃₄	34.27	N	14.01	(OCOCH ₃) ₅	295.22
C ₂₈	336.28	H ₃₅	35.28	N ₂	28.02	(OCOCH ₃) ₆	354.26
C ₂₉	348.29	H ₃₆	36.29	N ₃	42.03	(OCOCH ₃) ₇	413.31
C ₃₀	360.30	H ₃₇	37.30	N ₄	56.04	(OCOCH ₃) ₈	472.35
C ₃₁	372.31	H ₃₈	38.30	N ₅	70.05	(OCOCH ₃) ₉	531.40
C ₃₂	384.32	H ₃₉	39.31	N ₆	84.06	(OCOCH ₃) ₁₀	590.44
C ₃₃	396.33	H ₄₀	40.32	S	32.06	(H ₂ O) _{0.5}	9.01
C ₃₄	408.34	H ₄₁	41.33	S ₂	64.12	H ₂ O	18.02
C ₃₅	420.35	H ₄₂	42.34	S ₃	96.18	(H ₂ O) _{1.5}	27.02
H	1.008	H ₄₃	43.34	S ₄	128.24	(H ₂ O) ₂	36.03
H ₂	2.016	H ₄₄	44.35	F	19.00	(H ₂ O) ₃	54.05
H ₃	3.024	H ₄₅	45.36	F ₂	38.00	(H ₂ O) ₄	72.06
H ₄	4.032	H ₄₆	46.37	F ₃	57.00	(H ₂ O) ₅	90.08
H ₅	5.040	H ₄₇	47.38	Cl	35.45	(H ₂ O) ₆	108.10
H ₆	6.048	H ₄₈	48.38	Cl ₂	70.91		
H ₇	7.056	H ₄₉	49.39	Cl ₃	106.37		

www.alfa.com

NORTH AMERICA
Tel: +1-800-343-0660 or
+1-978-521-6300
Fax: +1-800-322-4757
Email: info@alfa.com

GERMANY
Tel: 00800 4566 4566 or
+49 721 84007 280
Fax: 00800 4577 4577 or
+49 721 84007 300
Email: Eurosales@alfa.com

UNITED KINGDOM
Tel: 0800-801812 or
+44 (0)1524-850506
Fax: +44 (0)1524-850608
Email: UKsales@alfa.com

FRANCE
Tel: 0800 03 51 47 or
+33 (0)3 8862 2690
Fax: 0800 10 20 67 or
+33 (0)3 8862 6864
Email: frventes@alfa.com

INDIA
Tel: +91 8008 812424 or
+91 8008 812525 or
+91 8008 812626
Fax: +91 8418 260060
Email: India@alfa.com

CHINA
Tel: +86 (010) 8567-8600
Fax: +86 (010) 8567-8601
Email: saleschina@alfa-asia.com

KOREA
Tel: +82-2-3140-6000
Fax: +82-2-3140-6002
Email: saleskorea@alfa-asia.com